

ARHGEF16 Polyclonal Antibody

Cat No: HR1AP6658

For research use only

Overview

Product Name	ARHGEF16 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	ARHGEF16
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARHGEF16; EPHEXIN4; NBR; Rho guanine nucleotide exchange factor 16; Ephexin-4
Background	Rho guanine nucleotide exchange factor 16(ARHGEF16) Homo sapiens Although the specific function of this protein is not known yet, it is thought to be involved in protein-protein and protein-lipid interactions. [provided by RefSeq, Jul 2008],