

ARHGAP9 Polyclonal Antibody

Cat No: HR1AP8275

For research use only

Overview

Product Name	ARHGAP9 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	90kDa
GeneID?Human?	ARHGAP9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARHGAP9; Rho GTPase-activating protein 9; Rho-type GTPase-activating protein 9
Background	Rho GTPase activating protein 9(ARHGAP9) Homo sapiens This gene encodes a member of the Rho-GAP family of GTPase activating proteins. The protein has substantial GAP activity towards several Rho-family GTPases in vitro, converting them to an inactive GDP-bound state. It is implicated in regulating adhesion of hematopoietic cells to the extracellular matrix. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008].