

ARHGAP12 Polyclonal Antibody

Cat No: HR1AP9147

For research use only

Overview

Product Name	ARHGAP12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	ARHGAP12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARHGAP12; Rho GTPase-activating protein 12; Rho-type GTPase-activating protein 12
Background	Rho GTPase activating protein 12(ARHGAP12) Homo sapiens This gene encodes a member of a large family of proteins that activate Rho-type guanosine triphosphate (GTP) metabolizing enzymes. The encoded protein may be involved in suppressing tumor formation by regulating cell invasion and adhesion. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012],