

RRP22 Polyclonal Antibody

Cat No: HR1AP5638

For research use only

Overview

Product Name	RRP22 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	RASL10A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RASL10A; RRP22; Ras-like protein family member 10A; Ras-like protein RRP22; Ras-related protein on chromosome 22
Background	similarity:Belongs to the small GTPase superfamily. Ras family.,