

# RASSF6 Polyclonal Antibody

Cat No: HR1AP4441

For research use only

## Overview

Product Name	RASSF6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	RASSF6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RASSF6; Ras association domain-containing protein 6
Background	<p>Ras association domain family member 6(RASSF6) Homo sapiens This gene encodes a member of the Ras-association domain family (RASSF). Members of this family form the core of a highly conserved tumor suppressor network, the Salvador-Warts-Hippo (SWH) pathway. The protein encoded by this gene is a Ras effector protein that induces apoptosis. A genomic region containing this gene has been linked to susceptibility to viral bronchiolitis. Alternative splicing results in multiple transcript variants and protein isoforms. [provided by RefSeq, Jul 2012].</p>