

# RASSF4 Polyclonal Antibody

Cat No: HR1AP4440

For research use only

## Overview

Product Name	RASSF4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	RASSF4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RASSF4; AD037; Ras association domain-containing protein 4
Background	Ras association domain family member 4(RASSF4) Homo sapiens The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described. [provided by RefSeq, Jul 2008],