

Rap 2C Polyclonal Antibody

Cat No: HR1AP4436

For research use only

Overview

Product Name	Rap 2C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	RAP2C
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAP2C; Ras-related protein Rap-2c
Background	RAP2C, member of RAS oncogene family(RAP2C) Homo sapiens The protein encoded by this gene is a member of the Ras-related protein subfamily of the Ras GTPase superfamily. Members of this family are small GTPases that act as molecular switches to regulate cellular proliferation, differentiation, and apoptosis. This protein has been reported to activate in vitro transcriptional activity of the serum response element. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2012].