

## Ran BP-6 Polyclonal Antibody

Cat No: HR1AP4434

For research use only

### Overview

Product Name	Ran BP-6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	125kDa
GeneID?Human?	RANBP6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RANBP6; Ran-binding protein 6; RanBP6
Background	function:May function in nuclear protein import as nuclear transport receptor.,similarity:Belongs to the importin beta family.,similarity:Contains 7 HEAT repeats.,