

## Rab 31 Polyclonal Antibody

Cat No: HR1AP5695

For research use only

### Overview

Product Name	Rab 31 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	RAB31
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB31; RAB22B; Ras-related protein Rab-31; Ras-related protein Rab-22B
Background	RAB31, member RAS oncogene family(RAB31) Homo sapiens Small GTP-binding proteins of the RAB family, such as RAB31, play essential roles in vesicle and granule targeting (Bao et al., 2002 [PubMed 11784320]).[supplied by OMIM, Jul 2009],