

## Rab 18 Polyclonal Antibody

Cat No: HR1AP6385

For research use only

### Overview

Product Name	Rab 18 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	RAB18
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB18; Ras-related protein Rab-18
Background	RAB18, member RAS oncogene family(RAB18) Homo sapiens The protein encoded by this gene is a member of a family of Ras-related small GTPases that regulate membrane trafficking in organelles and transport vesicles. Knockdown studies in zebrafish suggest that this protein may have a role in eye and brain development. Mutations in this gene are associated with Warburg micro syndrome type 3. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jan 2012],