

ARHGEF19 Polyclonal Antibody

Cat No: HR1AP2809

For research use only

Overview

Product Name	ARHGEF19 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	90kDa
GeneID?Human?	ARHGEF19
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARHGEF19; Rho guanine nucleotide exchange factor 19; Ephexin-2
Background	Rho guanine nucleotide exchange factor 19(ARHGEF19) Homo sapiens Guanine nucleotide exchange factors (GEFs) such as ARHGEF19 accelerate the GTPase activity of Rho GTPases (see RHOA, MIM 165390).[supplied by OMIM, Dec 2008].