

# ARHGEF10 Polyclonal Antibody

Cat No: HR1AP2808

For research use only

## Overview

Product Name	ARHGEF10 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	152kDa
GeneID?Human?	ARHGEF10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARHGEF10; KIAA0294; Rho guanine nucleotide exchange factor 10
Background	Rho guanine nucleotide exchange factor 10(ARHGEF10) Homo sapiens This gene encodes a Rho guanine nucleotide exchange factor (GEF). Rho GEFs regulate the activity of small Rho GTPases by stimulating the exchange of guanine diphosphate (GDP) for guanine triphosphate (GTP) and may play a role in neural morphogenesis. Mutations in this gene are associated with slowed nerve conduction velocity (SNCV). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015],