

# Rab L2A Polyclonal Antibody

Cat No: HR1AP4405

For research use only

## Overview

Product Name	Rab L2A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30-34kDa
GeneID?Human?	RABL2A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RABL2A; Rab-like protein 2A
Background	RAB, member of RAS oncogene family-like 2A(RABL2A) Homo sapiens This gene is a member of the RAB gene family which belongs to the RAS GTPase superfamily. The proteins in the family of RAS-related signaling molecules are small GTP-binding proteins that play important roles in the regulation of exocytotic and endocytotic pathways. This gene maps to the site of an ancestral telomere fusion event and may be a subtelomeric gene. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015].