

ARHGAP18 Polyclonal Antibody

Cat No: HR1AP2807

For research use only

Overview

Product Name	ARHGAP18 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	ARHGAP18
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ARHGAP18; Rho GTPase-activating protein 18; MacGAP; Rho-type GTPase-activating protein 18
Background	Rho GTPase activating protein 18(ARHGAP18) Homo sapiens ARHGAP18 belongs to a family of Rho (see MIM 165390) GTPase-activating proteins that modulate cell signaling (Potkin et al., 2009 [PubMed 19065146]).[supplied by OMIM, Apr 2010],