

# Rab 33A Polyclonal Antibody

Cat No: HR1AP9138

For research use only

## Overview

Product Name	Rab 33A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	RAB33A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB33A; RABS10; Ras-related protein Rab-33A; Small GTP-binding protein S10
Background	RAB33A, member RAS oncogene family(RAB33A) Homo sapiens The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. It is GTP-binding protein and may be involved in vesicle transport. [provided by RefSeq, Jul 2008],