

# Rab 40B Polyclonal Antibody

Cat No: HR1AP5686

For research use only

## Overview

Product Name	Rab 40B Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	RAB40B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB40B; SEC4L; Ras-related protein Rab-40B; SOCS box-containing protein RAR; Protein Rar
Background	RAB40B, member RAS oncogene family(RAB40B) Homo sapiens The protein encoded by this gene has similarity to a yeast protein which suggests a role of the gene product in regulating secretory vesicles. [provided by RefSeq, Jul 2008],