

Rab 7L1 Polyclonal Antibody

Cat No: HR1AP4404

For research use only

Overview

Product Name	Rab 7L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	RAB7L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB7L1; Ras-related protein Rab-7L1; Rab-7-like protein 1
Background	similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Ubiquitous.,