

G3BP2 Polyclonal Antibody

Cat No: HR1AP3496

For research use only

Overview

Product Name	G3BP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	G3BP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	G3BP2; KIAA0660; Ras GTPase-activating protein-binding protein 2; G3BP-2; GAP SH3 domain-binding protein 2
Background	function:Probable scaffold protein that may be involved in mRNA transport .,PTM:Arg-457 and Arg-468 are dimethylated, probably to asymmetric dimethylarginine.,similarity:Contains 1 NTF2 domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,subunit:Binds to the N-terminal domain of I-kappa-B-alpha.,