

Sos 2 Polyclonal Antibody

Cat No: HR1AP8364

For research use only

Overview

Product Name	Sos 2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	150kDa
GeneID?Human?	SOS2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SOS2; Son of sevenless homolog 2; SOS-2
Background	SOS Ras/Rho guanine nucleotide exchange factor 2(SOS2) Homo sapiens This gene encodes a regulatory protein that is involved in the positive regulation of ras proteins. Mutations in this gene are associated with Noonan Syndrome-9. [provided by RefSeq, Jul 2016],