

GBP5 Polyclonal Antibody

Cat No: HR1AP5767

For research use only

Overview

Product Name	GBP5 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	66kDa
GeneID?Human?	GBP5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GBP5; Guanylate-binding protein 5; GBP-TA antigen; GTP-binding protein 5; GBP-5; Guanine nucleotide-binding protein 5
Background	similarity:Belongs to the GBP family.,