

## Rad GTPase Polyclonal Antibody

Cat No: HR1AP4413

For research use only

### Overview

Product Name	Rad GTPase Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	RRAD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RRAD; RAD; GTP-binding protein RAD; RAD1; Ras associated with diabetes
Background	similarity:Belongs to the small GTPase superfamily. RGK family.,subunit:Interacts with calmodulin preferentially to the inactive, GDP-form. Binds CAMKII which is capable of phosphorylating RAD in vitro.,tissue specificity:Skeletal and cardiac muscle, and lung. Lesser amounts in placenta and kidney. Also detected in adipose tissue. Overexpressed in muscle of type II diabetic humans.,