

GBP4 Polyclonal Antibody

Cat No: HR1AP5766

For research use only

Overview

Product Name	GBP4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	GBP4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GBP4; Guanylate-binding protein 4; GTP-binding protein 4; GBP-4; Guanine nucleotide-binding protein 4
Background	guanylate binding protein 4(GBP4) Homo sapiens Guanylate-binding proteins, such as GBP4, are induced by interferon and hydrolyze GTP to both GDP and GMP (Vestal, 2005 [PubMed 16108726]).[supplied by OMIM, Dec 2008],