

RHOBTB3 Polyclonal Antibody

Cat No: HR1AP6369

For research use only

Overview

Product Name	RHOBTB3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	RHOBTB3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RHOBTB3; KIAA0878; Rho-related BTB domain-containing protein 3
Background	Rho related BTB domain containing 3(RHOBTB3) Homo sapiens RHOBTB3 is a member of the evolutionarily conserved RHOBTB subfamily of Rho GTPases. For background information on RHOBTBs, see RHOBTB1 (MIM 607351).[supplied by OMIM, Apr 2004],