

# CCP2 Polyclonal Antibody

Cat No: HR1AP3000

For research use only

## Overview

|                       |  |
|-----------------------|--|
| Product Name          | CCP2 Polyclonal Antibody   |
| Source                | Rabbit   |
| Applications          | WB,IHC-p,ELISA   |
| Species Reactivity    | Human,Mouse  |
| Recommended Dilutions |  |
| Immunogen             |  |
| Species               | Rabbit   |
| Storage               | -20°C/1 year   |
| Isotype               |  |
| Clonality             |  |
| Concentration         | 1 mg/ml  |
| Observed band         | 104kDa   |
| GeneID?Human?         | AGBL2  |
| Human Swiss-Prot No.  |  |
| Cellular localization |  |
| Alternative Names     | AGBL2; CCP2; Cytosolic carboxypeptidase 2; ATP/GTP-binding protein-like 2  |
| Background            | cofactor: Binds 1 zinc ion per subunit.,function: May play a role in the processing of tubulin.,sequence caution: Translated as Lys.,similarity: Belongs to the peptidase M14 family., |