

## Vav2 (phospho Tyr142) Polyclonal Antibody

Cat No: HR1AP8600

For research use only

### Overview

Product Name	Vav2 (phospho Tyr142) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	101kDa
GeneID?Human?	VAV2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	VAV2; Guanine nucleotide exchange factor VAV2; VAV-2
Background	vav guanine nucleotide exchange factor 2(VAV2) Homo sapiens VAV2 is the second member of the VAV guanine nucleotide exchange factor family of oncogenes. Unlike VAV1, which is expressed exclusively in hematopoietic cells, VAV2 transcripts were found in most tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2008],