

IQGAP3 Polyclonal Antibody

Cat No: HR1AP3751

For research use only

Overview

Product Name	IQGAP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	180kDa
GeneID?Human?	IQGAP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IQGAP3; Ras GTPase-activating-like protein IQGAP3
Background	similarity:Contains 1 CH (calponin-homology) domain.,similarity:Contains 1 Ras-GAP domain.,similarity:Contains 4 IQ domains.,