

EF-2 (phospho Thr56) Polyclonal Antibody

Cat No: HR1AP6155

For research use only

Overview

Product Name	EF-2 (phospho Thr56) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	EEF2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EEF2; EF2; Elongation factor 2; EF-2
Background	eukaryotic translation elongation factor 2(EEF2) Homo sapiens This gene encodes a member of the GTP-binding translation elongation factor family. This protein is an essential factor for protein synthesis. It promotes the GTP-dependent translocation of the nascent protein chain from the A-site to the P-site of the ribosome. This protein is completely inactivated by EF-2 kinase phosphorylation. [provided by RefSeq, Jul 2008],