

Rab 26 Polyclonal Antibody

Cat No: HR1AP4916

For research use only

Overview

Product Name	Rab 26 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	RAB26
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB26; Ras-related protein Rab-26
Background	RAB26, member RAS oncogene family(RAB26) Homo sapiens Members of the RAB protein family, including RAB26, are important regulators of vesicular fusion and trafficking. The RAB family of small G proteins regulates intercellular vesicle trafficking, including exocytosis, endocytosis, and recycling (summary by Seki et al., 2000 [PubMed 11043516]).[supplied by OMIM, Nov 2010],