

HSC70 Rabbit Polyclonal Antibody

Cat No: HR1AP9669

For research use only

Overview

Product Name	HSC70 Rabbit Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70-72kDa
GeneID?Human?	HSPA8 HSC70 HSP73 HSPA10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)
Background	<p>heat shock protein family A (Hsp70) member 8(HSPA8) Homo sapiens This gene encodes a member of the heat shock protein 70 family, which contains both heat-inducible and constitutively expressed members. This protein belongs to the latter group, which are also referred to as heat-shock cognate proteins. It functions as a chaperone, and binds to nascent polypeptides to facilitate correct folding. It also functions as an ATPase in the disassembly of clathrin-coated vesicles during transport of membrane components through the cell. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011],</p>