

GluR-?1 Polyclonal Antibody

Cat No: HR1AP6750

For research use only

Overview

Product Name	GluR-?1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	GRID1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GRID1; KIAA1220; Glutamate receptor delta-1 subunit; GluR delta-1 subunit
Background	glutamate ionotropic receptor delta type subunit 1(GRID1) Homo sapiens This gene encodes a subunit of glutamate receptor channels. These channels mediate most of the fast excitatory synaptic transmission in the central nervous system and play key roles in synaptic plasticity.[provided by RefSeq, Jan 2009],