

CLIC6 Polyclonal Antibody

Cat No: HR1AP7747

For research use only

Overview

Product Name	CLIC6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	CLIC6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CLIC6; CLIC1L; Chloride intracellular channel protein 6; Parchorin
Background	chloride intracellular channel 6(CLIC6) Homo sapiens This gene encodes a member of the chloride intracellular channel family of proteins. The gene is part of a large triplicated region found on chromosomes 1, 6, and 21. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2015],