

NHE-6 Polyclonal Antibody

Cat No: HR1AP5593

For research use only

Overview

Product Name	NHE-6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	SLC9A6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC9A6; KIAA0267; NHE6; Sodium/hydrogen exchanger 6; Na(+)/H(+) exchanger 6; NHE-6; Solute carrier family 9 member 6
Background	solute carrier family 9 member A6(SLC9A6) Homo sapiens This gene encodes a sodium-hydrogen exchanger that is amember of the solute carrier family 9. The encoded protein localizes to early and recycling endosomes and may be involved in regulating endosomal pH and volume. Defects in this gene are associated with X-linked syndromic mental retardation, Christianson type. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Apr 2010].