

Na⁺ CP type II? Polyclonal Antibody

Cat No: HR1AP4000

For research use only

Overview

Product Name	Na ⁺ CP type II? Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	SCN2B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SCN2B; Sodium channel subunit beta-2
Background	sodium voltage-gated channel beta subunit 2(SCN2B) Homo sapiens The protein encoded by this gene is the beta 2 subunit of the type II voltage-gated sodium channel. The encoded protein is involved in cell-cell adhesion and cell migration. Defects in this gene can be a cause of Brugada Syndrome, atrial fibrillation, or sudden infant death syndrome. [provided by RefSeq, Jul 2015],