

FAM80B Polyclonal Antibody

Cat No: HR1AP3435

For research use only

Overview

Product Name	FAM80B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	RIMKLB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RIMKLB; FAM80B; KIAA1238; Beta-citryl-glutamate synthase B; N-acetyl-aspartyl-glutamate synthetase B; NAAG synthetase B; NAAGS; Ribosomal protein S6 modification-like protein B
Background	cofactor: Binds 2 manganese ions per subunit.,similarity: Belongs to the rimK family.,similarity: Contains 1 ATP-grasp domain.,