

FAM80A Polyclonal Antibody

Cat No: HR1AP3434

For research use only

Overview

Product Name	FAM80A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	RIMKLA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RIMKLA; FAM80A; N-acetylaspartyl-glutamate synthetase A; NAAG synthetase A; NAAGS; Ribosomal protein S6 modification-like protein A
Background	cofactor: Binds 2 manganese ions per subunit.,similarity: Belongs to the rimK family.,similarity: Contains 1 ATP-grasp domain.,