

Sodium Potassium ATPase alpha-1 (Phospho-Tyr260) Antibody

Cat No: HR1AP9727

For research use only

Overview

Product Name	Sodium Potassium ATPase alpha-1 (Phospho-Tyr260) Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	115kDa
GeneID?Human?	ATP1A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Sodium/potassium-transporting ATPase subunit alpha-1 (Na ⁺ /K ⁺ ATPase alpha-1 subunit) (EC 3.6.3.9) (Sodium pump subunit alpha-1)
Background	ATPase Na ⁺ /K ⁺ transporting subunit alpha 1(ATP1A1) Homo sapiens The protein encoded by this gene belongs to the family of P-type cation transport ATPases, and to the subfamily of Na ⁺ /K ⁺ -ATPases. Na ⁺ /K ⁺ -ATPase is an integral membrane protein responsible for establishing and maintaining the electrochemical gradients of Na and K ions across the plasma membrane. These gradients are essential for osmoregulation, for sodium-coupled transport of a variety of organic and inorganic molecules, and for electrical excitability of nerve and muscle. This enzyme is composed of two subunits, a large catalytic subunit (alpha) and a smaller glycoprotein subunit (beta). The catalytic subunit of Na ⁺ /K ⁺ -ATPase is encoded by multiple genes. This gene encodes an alpha 1 subunit. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2009],