

CNG-2 Polyclonal Antibody

Cat No: HR1AP3124

For research use only

Overview

Product Name	CNG-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	CNGA2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CNGA2; CNCA; CNCA1; CNGG2; Cyclic nucleotide-gated olfactory channel; Cyclic nucleotide-gated cation channel 2; Cyclic nucleotide-gated channel alpha-2; CNG channel alpha-2; CNG-2; CNG2
Background	cyclic nucleotide gated channel alpha 2(CNGA2) Homo sapiens The protein encoded by this gene represents the alpha subunit of a cyclic nucleotide-gated olfactory channel. The encoded protein contains a carboxy-terminal leucine zipper that mediates channel formation. [provided by RefSeq, Jan 2010],