

AQP5 Polyclonal Antibody

Cat No: HR1AP2801

For research use only

Overview

Product Name	AQP5 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	AQP5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AQP5; Aquaporin-5; AQP-5
Background	aquaporin 5(AQP5) Homo sapiens Aquaporin 5 (AQP5) is a water channel protein. Aquaporins are a family of small integral membrane proteins related to the major intrinsic protein (MIP or AQP0). Aquaporin 5 plays a role in the generation of saliva, tears and pulmonary secretions. AQP0, AQP2, AQP5, and AQP6 are closely related and all map to 12q13. [provided by RefSeq, Jul 2008],