

CACNB3 Polyclonal Antibody

Cat No: HR1AP9331

For research use only

Overview

Product Name	CACNB3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	50kDa
GeneID?Human?	CACNB3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CACNB3; CACNLB3; Voltage-dependent L-type calcium channel subunit beta-3; CAB3; Calcium channel voltage-dependent subunit beta 3
Background	calcium voltage-gated channel auxiliary subunit beta 3(CACNB3) Homo sapiens This gene encodes a regulatory beta subunit of the voltage-dependent calcium channel. Beta subunits are composed of five domains, which contribute to the regulation of surface expression and gating of calcium channels and may also play a role in the regulation of transcription factors and calcium transport. [provided by RefSeq, Oct 2011],