

GABAA R?1 Polyclonal Antibody

Cat No: HR1AP6545

For research use only

Overview

Product Name	GABAA R?1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	GABRG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GABRG1; Gamma-aminobutyric acid receptor subunit gamma-1; GABA(A) receptor subunit gamma-1
Background	<p>gamma-aminobutyric acid type A receptor gamma1 subunit(GABRG1) Homo sapiens The protein encoded by this gene belongs to the ligand-gated ionic channel family. It is an integral membrane protein and plays an important role in inhibiting neurotransmission by binding to the benzodiazepine receptor and opening an integral chloride channel. This gene is clustered with three other family members on chromosome 4. [provided by RefSeq, Jul 2008],</p>