

AChR α 1 Polyclonal Antibody

Cat No: HR1AP5751

For research use only

Overview

Product Name	AChR α 1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CHRNA1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CHRNA1; ACHRA; CHNRA; Acetylcholine receptor subunit alpha
Background	cholinergic receptor nicotinic alpha 1 subunit(CHRNA1) Homo sapiens The muscle acetylcholine receptor consists of 5 subunits of 4 different types: 2 alpha subunits and 1 each of the beta, gamma, and delta subunits. This gene encodes an alpha subunit that plays a role in acetylcholine binding/channel gating. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Nov 2012].