

# CLCC1 Polyclonal Antibody

Cat No: HR1AP3105

For research use only

## Overview

|                       |   |
|-----------------------|---|
| Product Name          | CLCC1 Polyclonal Antibody   |
| Source                | Rabbit  |
| Applications          | WB,IF,ELISA   |
| Species Reactivity    | Human   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 62kDa   |
| GeneID?Human?         | CLCC1   |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | CLCC1; KIAA0761; MCLC; Chloride channel CLIC-like protein 1; Mid-1-related chloride channel protein 1     |
| Background            | function:Seems to act as a chloride ion channel.,similarity:Belongs to the chloride channel MCLC family., |