

# LC3B Rabbit pAb

Cat No: HR1APEA174

For research use only

## Overview

Product Name	LC3B Rabbit pAb
Source	Rabbit
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:200-500
Immunogen	
Species	Rabbit
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1mg/ml
Observed band	14,16kDa
GeneID?Human?	81631
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ATG8F antibody, MAP1 light chain 3 like protein 2 antibody, MAP1A/MAP1B LC3 B antibody
Background	Autophagy marker Light Chain 3 (LC3) was originally identified as a subunit of microtubule-associated proteins 1A and 1B (termed MAP1LC3).