

Brp44 Polyclonal Antibody

Cat No: HR1AP6556

For research use only

Overview

Product Name	Brp44 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	MPC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MPC2; BRP44; Mitochondrial pyruvate carrier 2; Brain protein 44
Background	similarity:Belongs to the UPF0041 (0-44) family.,