



KV?1 Rabbit pAb

Cat No: HR1APEA307

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | KV?1 Rabbit pAb |
| Source | Rabbit |
| Applications | IHC |
| Species Reactivity | Human,Rat,Mouse |
| Recommended Dilutions | IHC 1:100-200 |
| Immunogen | |
| Species | Rabbit |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Concentration | 1mg/ml |
| Observed band | 50kDa |
| GeneID?Human? | 7881 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | hKvb3,KCNA1B,KV beta 1,Voltage-gated potassium channel beta 1 subunit |
| Background | N/A |