



Kv1.3 Rabbit pAb

Cat No: HR1APEA308

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | Kv1.3 Rabbit pAb |
| Source | Rabbit |
| Applications | WB |
| Species Reactivity | Human,Rat,Mouse |
| Recommended Dilutions | WB 1:1,000-2,000 |
| Immunogen | |
| Species | Rabbit |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Concentration | 1mg/ml |
| Observed band | 64kDa |
| GeneID?Human? | 3738 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | HGK5,KCNA3,MK3,Potassium channel 3, voltage gated potassium channel protein kv1.3 |
| Background | N/A |