

ACAD-11 Polyclonal Antibody

Cat No: HR1AP8846

For research use only

Overview

Product Name	ACAD-11 Polyclonal Antibody
Source	Rabbit
Applications	WB, IHC-p, ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	87kDa
GeneID?Human?	ACAD11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ACAD11; Acyl-CoA dehydrogenase family member 11; ACAD-11
Background	acyl-CoA dehydrogenase family member 11(ACAD11) Homo sapiens This gene encodes an acyl-CoA dehydrogenase enzyme with a preference for carbon chain lengths between 20 and 26. Naturally occurring read-through transcription occurs between the upstream gene NPHP3 (nephronophthisis 3 (adolescent)) and this gene. [provided by RefSeq, Aug 2015],