

ABHD9 Polyclonal Antibody

Cat No: HR1AP8709

For research use only

Overview

Product Name	ABHD9 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	EPHX3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EPHX3; ABHD9; Epoxide hydrolase 3; Abhydrolase domain-containing protein 9
Background	EPHX3 (Epoxide Hydrolase 3) is a Protein Coding gene. GO annotations related to this gene include hydrolase activity. An important paralog of this gene is EPHX4.