

ABHD8 Polyclonal Antibody

Cat No: HR1AP8692

For research use only

Overview

Product Name	ABHD8 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	ABHD8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ABHD8; Abhydrolase domain-containing protein 8
Background	abhydrolase domain containing 8(ABHD8) Homo sapiens This gene is upstream of, and in a head-to-head orientation with the gene for the mitochondrial ribosomal protein L34. The predicted protein contains alpha/beta hydrolase fold and secretory lipase domains. [provided by RefSeq, Jul 2008],