

ABHD14A Polyclonal Antibody

Cat No: HR1AP6553

For research use only

Overview

Product Name	ABHD14A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	ABHD14A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ABHD14A; Abhydrolase domain-containing protein 14A
Background	caution:It is uncertain whether Met-1 or Met-30 is the initiator.,function:Possible role in granule neuron development.,similarity:Belongs to the AB hydrolase superfamily. ABHD14 family.,